



Raz-AP

Attendant Propel

Australian Order Forms
Effective April 1, 2024

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AP: Attendant Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 185kg.**

1. Standard Raz-AP Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]

Z10018	Raz-AP w/ 18"W base frame – Standard base frame width, will fit over most toilets	ADD
--------	---	-----

2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.3]

Z10016	Raz-AP w/ 16"W frame width ¹ (base chair code and custom upcharge code below will be added)	ADD
Z1SW16	Upcharge for custom 16"W frame – will not fit over most toilets	ADD
Z10020	Raz-AP w/ 20"W frame width ² (base chair code and custom upcharge code below will be added)	ADD
Z1SW20	Upcharge for custom 20"W frame	ADD
Z10022	Raz-AP w/ 22"W frame width (base chair code and custom upcharge code below will be added)	ADD
Z1SW22	Upcharge for custom 22"W frame	ADD
Custom Frame Configurations		
ZLNGF3	Upcharge for custom 3" Longer Frame rec. for seat depths 20"/21", req for seat depths 21"-24".	ADD
Z1LOWF4	Upcharge for custom Low Frame (16"-18" STF, 4" lower than standard)	ADD
Z1APRF3	Upcharge for Recliner Frame (3" longer at the rear for added stability)	ADD
¹ 18"W Visco-Foam Interface seats can be modified to fit.		
² 18"W seats can be selected but must also select Seat Support Kit ZSSKAS [Sec.3]		

3. Adjusta-Back Frames – Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 4

Includes height- and angle-adjustable arms, tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, on which the upholstery cannot be raised.

ZAB __ 15*	Adjusta-Back ¹ – 15"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZAB __ 17*	Adjusta-Back ¹ – 17"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZAB __ 20*	Adjusta-Back ¹ – 20"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZAB __ 24*	Adjusta-Back ¹ – 24"H	* insert back frame width to complete the code (16" / 18" / 20" / 22" / 24")	ADD


Reclining Backs – 30° (must select Reclining Frame option Z1APRF3 in Sec 2) Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLAS / Z2PSAS) are included with Reclining Backs.

ZRBA __ 20*	Reclining Adjusta-Back with 20"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	ADD
ZRBA __ 24*	Reclining Adjusta-Back with 24"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	ADD
ZDTL	Dual Triggers – both on left side	N/C
ZDTR	Dual Triggers – both on right side	N/C

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation.

4. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.3] (configurations below are not available for recliner backs)

NOTE: For back frame width variances on recliner frames, please request custom quotation.

ZNB __ __ *	Custom – Back Frame 2" Narrower than chair frame	ADD	Chair Frame Width	16"	18"	20"	22"	
* Enter chair and back frame widths to complete code. eg. ZNB2018 specifies a 20"W chair frame with an 18"W back frame			Back Frame Width	14" ¹	16" ²	18" ³	20" ⁴	

¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit Z4SSKAS.

⁴ Custom Seat recommended.

5. Gel Back Support - Layer Gel Pressure reduction back upholstery for Raz Backrest. **New**

ZGB __ 15*	Gel Back – 15"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZGB __ 17*	Gel Back – 17"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZGB __ 20*	Gel Back – 20"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD
ZGB __ 24*	Gel Back – 24"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	ADD

6. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

ZSB14	14"W Symphony Back ¹ – recessed solid back support	ADD
ZSB16	16"W Symphony Back ¹ – recessed solid back support	ADD
ZSB18	18"W Symphony Back ¹ – recessed solid back support	ADD
ZSB20	20"W Symphony Back ¹ – recessed solid back support	ADD
ZSB22	22"W Symphony Back ¹ – recessed solid back support	ADD
ZHB	Harmony Back ^{1,2} (contoured solid back – 18"W x 13"H) Fits 18"W back frames only Fits 17" / 20" back heights only	ADD

¹ Does not fit 15"H backs

² Not recommended for 24"H backs or reclining backs.

7. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16"D and 19"D, are now designated 17"D and 20"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
ZMS17FRO	18"W x 17"D – Front Access	ADD	ZVFI17BRI	18"W x 17"D – Front Bridged	ADD
ZMS17BRI	18"W x 17"D – Front Bridged	ADD	ZVFI17FRO	18"W x 17"D – Front Access	ADD
ZMS20FRO	18"W x 20"D – Front Access	ADD	ZVFI17LTO	18"W x 17"D – Left Access	ADD
ZMS20BRI	18"W x 20"D – Front Bridged	ADD	ZVFI17RTO	18"W x 17"D – Right Access	ADD
Specialty Seats			ZVFI17RRO	18"W x 17"D – Rear Access	ADD
ZSOM	Shower-Only Seat – Molded 18"W x 16"-21"D	ADD	ZVFI20BRI	18"W x 20"D – Front Bridged	ADD
ZEJRM	E&J Replica Seat – Molded – 18"W x 17"D	ADD	ZVFI20FRO	18"W x 20"D – Front Access	ADD
ZEJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	ADD	ZVFI20LTO	18"W x 20"D – Left Access	ADD
			ZVFI20RTO	18"W x 20"D – Right Access	ADD
			ZVFI20RRO	18"W x 20"D – Rear Access	ADD
Seat Options					
ZSSKAS	Seat Support Kit for Raz-AP - These are used to support the sides of seats that are narrower than seat frames on which they are mounted. The reason for a narrower seat on a wider frame is typically to select a Molded Seat and avoid the upcharge of a custom seat. A typical example would be when a 20"W chair is needed but an 18"W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.				ADD
ZSS05	Spacers for Front Seat Brackets – 1/2" (pr) – for adding dump angle (based on 18"D seat)				ADD
ZSS10	Spacers for Front Seat Brackets – 1" (pr) – for adding dump angle (based on 18"D seat)				ADD

8. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

ZCSEAT	Custom seat Select if both width and depth are <19"	ADD
ZCSEATHD	Custom seat Select if either width or depth are between 20"<24"	ADD
ZCSEATHDX	Heavy Duty Custom seat Select if width and/or depth is ≥ 25"	ADD

Custom Gel Sheet **New**

GRSPTS2428	3/4" Pillow Top Gel Scrimmed upgrade to Custom Seat. (Replaces Visco Foam with gel. Up to dimensions of 20" x 24")	ADD
------------	---	-----

MUST FILL OUT PARAGON CUSTOM SEAT FORM ON BACK PAGE

- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

9. Commode Set up - **These set up configurations are for dealer use only. Commodes will not be supplied out of the box with the below measurements.**

Set armrest height to _ _ _ _ Measured from underside of the seat to the top of armrest.

Seat toilet clearance to _ _ _ _ This is the minimum clearance between the rear bar of the commode and the ground. Excludes Zimpankit)

Footrest length _ _ _ _ Measured from underside of the seat to the top of the footplate.

10. Arm Supports and Accessories

ZARML	Flip-up Arm Supports (pr) Long w/ 13" pad	STD	ZPAML	Pivoting Arm Mount – Left (for ZFPLL/ZTRL)	ADD
ZARMS	Flip-up Arm Supports ¹ (pr) Short w/ 10" pad	N/C	ZPAMR	Pivoting Arm Mount – Right (for ZFPLR / ZTRR)	ADD
ZALK	Arm Support Locks (pr)	ADD	ZASK05	Arm Support Spacer Kit – 1/2" per side for all Adjusta-Backs (pr)	ADD
Z2PLAS	2-Point Arm Supports ² (pr) – Long with 13"L pads (std with reclining backs)	ADD	ZASK10	Arm Support Spacer Kit – 1" per side for all Adjusta-Backs (pr)	ADD
Z2PSAS	2-Point Arm Supports ² (pr) – Short with 10"L pads (std with reclining backs)	ADD	ZPHGL	Pivoting Hand Grips (pr) for 13"L arms ⁴	ADD
ZFAPLL	Long Molded Flat Armpad ³ – Left (14"L x 3.25"W – fits ZARM13)	ADD	ZPHGS	Pivoting Hand Grips (pr) for 10"L arms ⁴	ADD
ZFAPLR	Long Molded Flat Armpad ³ – Right (14"L x 3.25"W – fits ZARM13)	ADD	ZAPSB_ _ *	Anterior Posture Support Bar ⁵ * Enter frame width to complete code (14" / 16" / 18" / 20")	ADD
ZFAPSL	Short Molded Flat Armpad ³ – Left (10"L x 3"W – fits ZARM10)	ADD	ZAPSB22	Anterior Posture Support Bar ⁵ for a 22"W Back Frame	ADD
ZFAPSR	Short Molded Flat Armpad ³ – Right (10"L x 3"W – fits ZARM10)	ADD			
ZTRL	Molded Arm Trough ³ – Left (18.5"L x 4.5"W)	ADD			
ZTRR	Molded Arm Trough ³ – Right (18.5"L x 4.5"W)	ADD			
¹ Not compatible with ZFAPL / ZFAPS / ZTR / ZPAM			⁴ ZPHGL / ZPHGS combined with 2 Point Arms requires a custom upcharge.		
² For chairs with same back width and frame width only			⁵ Not compatible with ZFAPL / ZFAPS / ZTR		
³ May require ZPAM / ZASK05 / ZASK10 to flip up, depending on medial / lateral position of the Flat Arm Pad or Arm Trough					

11. Foot and Leg Supports

ZFS	Foot Supports (pr) – 16"-20.5"	STD	ZELRL	Elevating Leg Support ³ – Left (includes ZADML)	ADD
ZFS16	Foot Supports for 16"W frames (pr) – (16"-20.5")	N/C	ZELRR	Elevating Leg Support ³ – Right (includes ZADMR)	ADD
ZFET2	Footrest Extension Tubes ¹ (pr) – 2" longer	ADD	ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 __	ADD
ZMFX8	MFX8 Foot Supports (pr) – 8"- 12"	ADD	ZRLSL	Universal Residual Limb Support – Left	ADD
ZMFX12	MFX12 Foot Supports (pr) – 12"-16.5"	ADD	ZRLSR	Universal Residual Limb Support – Right	ADD
ZADM	Adj Angle/Depth Molded Footplates ² (pr)	ADD	ZFBAS	Flip-Back Foot Support (1-piece) 16.5"-20"	ADD
ZFET2A	Footrest Extension Tubes (pr) for ZADML – 2" longer	ADD			
ZADABS__*	Angle/Depth Adj ABS Footplates (pr) * Enter frame width to complete code	ADD			
ZVS__*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (16" / 18")	ADD			
ZVS__*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (20" / 22" / 24")	ADD			

¹ Not compatible with ZMFX12/ZADM/ZVS/ZFBT

² Not compatible with Z30016 or Z7MFX8

³ ZELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 3" longer seat frames (Z3LNGSF3)

12. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	ADD	ZNC	Neoprene Cover for Standard Footplate (pr)	ADD
ZCS__*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20" / 22" / 24")	ADD	ZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates	ADD
ZCSIC__*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	ADD	ZNC16	Neoprene Covers for Narrow Molded Footplates – 16"W frame (pr)	ADD
ZCSIC__*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (22" / 24")	ADD	ZNCABS__*	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable Footplates	ADD
ZCPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths (16" / 18" / 20")	ADD	ZNCVS__*	Neoprene Cover for V-Style Foot Support	ADD
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	ADD	ZNCFB	Neoprene Cover for Flip-Back Foot Support	ADD
ZHS2224	H-Strap – for 22" / 24"W frames – hook and loop attachment	ADD	ZLMOR	Lateral / Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 __	ADD
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	ADD	ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 __	ADD

¹ Not compatible with ZMFX12 / ZFBT / ZVS

13. Gel Accessories New

GHPSTD	Geo - Matrix Hanger Pads Standard Size Size 6' x 8' (Pair)	ADD
GHPLRG	Geo - Matrix Hanger Pads Large Size Size 8.5' x 10.5' (Pair)	ADD
GFPSTD	Geo - Matrix Foot Plate Pads Standard Size "Cascading" 6' x 8' (Pair)	ADD
GFPLRG	Geo - Matrix Foot Plate Pads Large Size "Cascading" 8' x 10' (Pair)	ADD

14. Positioning Pads - Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

Head Supports		
ZHSAS1	Head Support Pad – Small (10" x 4")	ADD
ZHSASA2	Head Support Pad – Large (12" x 6")	ADD
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	ADD
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate. (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	ADD
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	ADD
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	ADD
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	ADD
Leg Adductors (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
ZLGT1L	Leg Adductor w/ Small Pad – Left	ADD
ZLGT1R	Leg Adductor w/ Small Pad – Right	ADD
ZLGT2L	Leg Adductor w/ Large Pad – Left	ADD
ZLGT2R	Leg Adductor w/ Large Pad – Right	ADD
Hip Adductors (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
ZHP1L	Hip Adductor w/ Small Pad – Left	ADD
ZHP1R	Hip Adductor w/ Small Pad – Right	ADD
ZHP2L	Hip Adductor w/ Large Pad – Left	ADD
ZHP2R	Hip Adductor w/ Large Pad – Right	ADD
Hip Adductors For Reclining Backs (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left	ADD
ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right	ADD
ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left	ADD
ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right	ADD
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors These plates allow for an additional 1.5" medial offset and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
ZLAEP	Lateral/Adductor Extension Plate (ea) Total QTY __	ADD
Medial Thigh Supports (Pommel) (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
ZMTST1	Medial Thigh Support w/ Small Pads	ADD
ZMTST2	Medial Thigh Support w/ Large Pads	ADD
Neoprene Covers for Positioning Pads		
☐ PC-NCP1	Neoprene cover for Positioning Pad -Small Total QTY __	ADD
☐ PC-NCP2	Neoprene cover for Positioning Pad - Large Total QTY __	ADD
☐ PC-GNCP1	Neoprene cover with Gel Insert for Positioning Pad - Small New Total QTY __	ADD
☐ PC-GNCP2	Neoprene cover with Gel Insert for Positioning Pad - Large New Total QTY __	ADD

15. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!


Pelvic Belts			Chest Belts		
3001	Kinetic Balance Positioning Products Belt - S New	ADD	3001	Kinetic Balance Positioning Products Belt - S New	ADD
3002	Kinetic Balance Positioning Products Belt - M New	ADD	3002	Kinetic Balance Positioning Products Belt - M New	ADD
3003	Kinetic Balance Positioning Products Belt - L New	ADD	3003	Kinetic Balance Positioning Products Belt - L New	ADD
ZPBS	Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	ADD	ZCBS	Chest Belt ² – Small (back widths 14"- 16", 1-piece w/ hook and loop closure)	ADD
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	ADD	ZCBM	Chest Belt ² – Medium (back widths 16"- 18", 1-piece w/ hook and loop closure)	ADD
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	ADD	ZCBL	Chest Belt ² – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	ADD
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22"- 26", 2-piece w/ side-release buckle)	ADD	ZCBXL	Chest Belt – XLarge ² (back frame widths 22"-26", 1-piece w/ hook and loop closure)	ADD
ZPBICS	Infection Control Pelvic Belt ¹ – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	ADD	ZCBICS	Infection Control Chest Belt ^{1,2} – Small (back widths 14"-16", 1-piece w/ side-release buckle)	ADD
ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	ADD	ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	ADD
ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	ADD	ZCBICL	Infection Control Chest Belt ^{1,2} – Large (back widths 18"-22", 1-piece w/ side-release buckle)	ADD
ZPBICXL	Infection Control Pelvic Belt ¹ – XLarge (seat frame widths 22"-26", 2-piece w/ side-release buckle)	ADD	ZCBICXL	Infection Control Chest Belt ^{1,2} – XLarge (back frame widths 22"-26", 1-piece w/ side-release buckle)	ADD

¹ Non-absorbent belts – recommended for infection-control environments

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM / ZCBL / ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

16. Kinetic Conversion Information **New**

Old Bodypoint Code & Size	Use	New Kinetic Size - See size Chart
ZPBBPS - Small	→	3001 - Small
ZPBBPM - Medium	→	3001 - Small
ZPBBPL - Large	→	3002 - Medium
ZPBBPXL - Extra Large	→	3002 - Medium
ZPBBPXLHD - Extra Large HD	→	3003 - Large (Include spacers and Brackets)
ZPBBPXLHDAS - Extra Large HD (AP/SP600)	→	3003 - Large (Include spacers and Brackets)



	A	B	C
S	330	110	38
M	430	110	38
L	510	110	38

All measurements in mm

17. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	ADD
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	ADD
ZSC17	Splash Curtain ¹ for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	ADD
ZSC20	Splash Curtain ¹ for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	ADD
ZSCFR	Splash Curtain for 17" / 20"D (at front of aperture only, not for Front Open Seats)	ADD
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	ADD
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	ADD
ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	ADD

¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.

18. Other

ZDCLA	Directional Caster Lock – fits on left rear caster (facilitates steering)	ADD	ZPSK1	Commode Pan Spacer Kit ² – 1" lower	ADD
ZATAS	Anti-Tippers ¹ (pr) Anti-Tippers are sized according to caster diameter ordered. Rear Front	ADD	ZPSK2	Commode Pan Spacer Kit ² – 2" lower	ADD
			ZPSK3	Commode Pan Spacer Kit ² – 3" lower	ADD
			ZPSK4	Commode Pan Spacer Kit ² – 4" lower	ADD
ZCUAP4	4" Caster Upgrade (4x)	ADD	ZIMPANKIT	Commode Pan and Seat attachments	ADD
ZCUAP6	6" Caster Upgrade (4x)	ADD	ZSCLAMPKIT	Seat Clamp Kit. (Only required if not ordering a seat from Paragon)	ADD

¹ Rear Anti-Tippers are included with Reclining Back Option [Sec.7]

² Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

Some assembly of Raz products is required.

19. Notes:

Specifications

Raz-AP (Attendant Propel)

Weight Capacity	185kg
Overall Width	Seat width plus 4.5" (with Standard Arm Support Pads)
Overall Length (foot supports in middle position)	36.5" (with rear caster trailing), 35" (to rear canes)
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat)	
- with pan brackets	15.25"-19.25" (1" increments)
- without pan brackets	16"-20" (1" increments)
Toilet Clearance Width	16" for standard 18"W seat frame. For other widths, seat frame width minus 2"
Arm Support Height - above molded seat (1/2" increments)	
- 15"H Adjusta-Back	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5" (1/2" increments)
- 20"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 20"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	20"-24"
Product Weight (with 18"W seat)	17kg
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)

Paragon - CUSTOM SEAT ASSESSMENT FORM



Client:

Reference #:

- ☐ ZCSEAT **Select if both width and depth are <19"**
- ☐ ZCSEATHD **Select if either width or depth are between 20"<24"**
- ☐ ZCSEATHDX **Select if width and/or depth is ≥ 25"**

Custom Gel Sheet **New**

- ☐ GRSPTS2428 3/4" Pillow Top Gel Scrimmed upgrade to Custom Seat.
(Replaces Visco Foam with gel. Up to dimensions of 20" x 24")


Chair Model:

- ☐ AT ☐ CAT ☐ ART
- ☐ ZUMAP/SP ☐ AP/SP
- ☐ AP/SP600 ☐ AT600




Serial #: _____


Section 1. Access Opening: **Must select one.**




☐ Bridged Front




☐ Front Open



☐ Left Open




☐ Right Open




☐ Rear Open


Section 2. Aperture Shape: **Must select one.**




☐ No Aperture





☐ Keyhole - Style



☐ CAT - Style



☐ EJR - Style



- Notes:
1. Keyhole Style Aperture Shown
2. All units to be in inches
3. Seat width can be up to 2" wider than seat frame
4. Call for other customization options
5. Reference to a "standard" seat refers to an 18" wide Seat dimension

A	Seat Width:	E	Aperture Width:
B	Seat Depth:	F	Aperture Depth:
C	Seat Frame Width:	G	Opening Width:
D	Back Frame Width:	H	Aperture Location:

Seat Frame Depth **(Must Select one):** ☐ Std ☐ +3" (Z3LNGSF3) ☐ Other (Specify)

Comments: