

Universal Positioning Bases


Traditionally, it has been difficult to meet clinical requirements for wheelchair users with pelvic and postural requirements when seated on pressure care cushions. The Universal Positioning Base solves these challenges by insertion under the air or gel cushion to configure the base system as per the user's clinical assessment to achieve a functional outcome whilst not compromising with pressure care properties.

ENHANCED POSITIONING OUTCOMES UNDERNEATH AIR AND GEL CUSHIONS




Can be **INSERTED INTO** the existing cushion cover where suitable depending on caregiver requirement.

ATTACHES to the wheelchair seat pan and remains attached to the wheelchair at all times.



POSITIONING FOR GEL CUSHIONS



POSITIONING FOR AIR CUSHIONS


Base Width (inches)		14	15	16	17	18	19	20
Base Length (inches)	14							
	15							
	16							
	17							
	18							
	19							
	20							



CAN BE CUT AND TRIMMED

- To order Base Options and Accessories replace **WW** in the code with Width and **LL** with Length


BASE OPTIONS



CUSHION SPACER KIT

Raises cushion for extra long leg lengths, or use as a base for attaching contouring support wedges.


1113-01-WWLL



ISCHIAL CONTOUR CUSHION SPACER KIT

Helps to prevent posterior pelvic sliding.

1113-02-WWLL



HIP FLEXION CONTOUR SPACER KIT

Helps to prevent posterior pelvic sliding. Allows for 10° Hip Flexion relief.

1113-03-WWLL

ACCESSORIES



CONTOUR KIT

Provides additional support for pelvis and thighs. Includes 2 lateral thigh contour wedges and a medial thigh contour wedge.

1113-05-00LL



PRESSURE REDUCING COVER

Includes neoprene outer cover and pressure-reducing top layer. Provides back-up pressure relief if air cushions deflate.

1113-07-WWLL